

Tailor your system



The package includes many useful accessories. If you want a solution completely tailored to your needs, you can also select from a wide range of original accessories. Choose Domino Classic and discover the easy way to interact with your friends, colleagues and classmates.

The package includes

- Transmitter unit
- Receiver unit
- Durable travel case
- Universal charger kit
- Neck straps & belt clips
- Earphones or Neckloop



Discover Domino Classic

When was the last time you heard music and speech the way it's meant to sound?

Domino Classic is a listening system that works together with your hearing aids to deliver every word in its natural essence – even in the toughest listening situation. It's the perfect choice in circumstances where your hearing aids just aren't enough, but it works just as well for people who just suffer the occasional hearing difficulty.

Listening to a distant speaker at a seminar, enjoying the company of friends at a busy restaurant or watching a TV show – Domino Classic delivers speech with incredible detail. Directly in your ears and down to the last syllable.

Packed with the latest digital technology it's still surprisingly easy to use.



Dedicated to your hearing

Bellman & Symfon of Sweden is dedicated to improve the quality of life for people with hearing loss. Our people and partners are devoted to this mission and we work closely with hearing healthcare professionals and leading experts to protect your hearing and make better hearing possible for everyone.

European office

Bellman & Symfon Europe AB
Södra Långebergsgatan 30
421 32 Västra Frölunda
Sweden
Phone +46 31 68 28 20
Fax +46 31 68 28 90
E-mail info@bellman.se


Asia Pacific office

Bellman & Symfon Asia Ltd.
1401-06, Main Tower, GD Int'l Bldg
339 Huanshi Dong Rd
Guangzhou, China, 510095
Phone +86 20 8318 2950
Fax +86 20 8333 5528
E-mail info.cn@bellman.com

© 2017, Bellman & Symfon AB, all rights reserved. Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB.

bellman.com

Domino Classic Personal listening system



Are these situations familiar?

-  **At home.** Is it hard to follow the conversation when you are talking to friends and relatives?
-  **At work.** Do you have trouble hearing what's being said at busy meetings or seminars?
-  **At school.** Do you find it difficult to hear what your teacher is saying in class?
-  **Watching TV.** Do you miss words when you are watching TV at normal volume?
-  **Enjoying music.** Do your favorite songs not sound as good as they used to?
-  **Using hearing aids.** Do you feel that your hearing aids don't fully help in certain situations?

Discover the **easy** way to hear better



High Definition Audio

Domino Classic uses the latest digital technology to stream audio in full stereo. The result is speech and music with outstanding clarity and definition.



Wireless Secure Encryption

Domino Classic uses 128-bit encryption, just like your computer at work. This means you can use your system even at confidential business meetings.



Travel Smart 2.4 GHz

Domino Classic uses the 2.4 GHz global frequency standard. This means you can use your Domino system anywhere in the world.



Easy Connect







Domino Classic is ready for action with the touch of a button, and you can just as easily connect it to your TV, computer, stereo system or MP3 player.



 **Bellman & Symfon®**
DESIGN FOR EARS

CL_001854001_EN

Discover the **easy** way to hear better




How Domino Classic works

The system consists of a transmitter and a receiver, both equipped with state-of-the-art microphones.

How it works. Sound is picked up by the microphones and converted into digital data. Speech is clarified and noise is eliminated. The speech is then streamed over an encrypted broadband link – from the transmitter to the receiver. The receiver amplifies the speech and sends it to your headphones, earphones or neck loop. From the receiver, you can easily control the functions of the two units.

More than words. The transmitter also incorporates a stereo input for streaming high quality audio from your TV, computer, MP3 player or stereo system.



Options for every need

<p>EARPHONES</p>	<p>HEADPHONES</p>	<p>CAR CHARGER</p>
<p>NECK LOOP</p> <p>NEW!</p>	<p>These original accessories are especially developed to maximize the performance of your Domino Classic system. Contact your nearest dealer or visit www.bellman.se/audio</p>	
<p>TIE CLIP MIC.</p>	<p>CABLE KIT, 16" (5m)</p>	<p>EXTRA RECEIVER</p>